



This PDF is generated from authoritative online content, and is provided for convenience only. This PDF cannot be used for legal purposes. For authoritative understanding of what is and is not supported, always use the online content. To copy code samples, always use the online content.

# Performance DNA Administrator Guide

Forget archived users

12/15/2025

## Forget archived users

This option is available on the **System Settings > Event Settings** tab. It allows you to set the frequency to run the forget archived users event. To enable it, select the **Enabled** check box for the **Forget archived users** item. The **Interval in Minutes** setting specifies how frequently you want to run the **Forget archived users** event.

Note that this event depends on the value of the **Number of days after which archived users are eligible for auto-anonymisation** field in the **General Settings** tab. The **archived users** who exceed the number of days specified in the aforementioned field will be forgotten automatically when the system runs this event.